



Reagent Kit for EZ1010 Copper Cu(II) Analyzer

de producto: **APPC1010KTO**
USD Precio: Contacto Hach Chile
Llamar para confirmar fecha de envío

Especificaciones

Contenido de la caja: APPC101001 Nitric Acid Solution, 0.25 M, 2 L (2x)
APPC230201 EZ2302, Color Indicator Solution, 2 L (2x)
259342 Standard Solution, Copper, 1000 mg/L Cu (NIST), 100 mL (1x)

Parámetro: Copper, dissolved

Contenido de la caja

APPC101001 Nitric Acid Solution, 0.25 M, 2 L (2x) APPC230201 EZ2302, Color Indicator Solution, 2 L (2x) 259342 Standard Solution, Copper, 1000 mg/L Cu (NIST), 100 mL (1x)